

BOOT AND GLOVES COLLECTION SYSTEM**Presentation**

Pair of nitrile gloves and pair of sterile plastic boots.

Sterilization method

Gamma irradiation.

Application

Ideal tool to avoid contamination of the swab used in the aseptic collection of samples for environmental monitoring.

How to use

After putting on the appropriate clothing and sanitizing your hands, open the kit packaging, remove the plastic boots and put them on, then remove the gloves and put them on.

Precautions and special care

Product intended for *in vitro* diagnostic use only.
Restricted for use by professionals. Do not inhale or ingest.
Do not use the product beyond the expiration date, with signs of contamination, or if it has changed color. In the presence of contamination, the product should be immediately discarded.
Do not use the product if the packaging is damaged or tampered with.

Storage

Store at room temperature in a dry place and protect from light.

Shelf-life

2 years from the date of manufacture.

Disposal of the product

After use, the product must be handled at the generating unit before environmentally appropriate final disposal, in accordance with official regulations.

Quality Guarantee

bioBoaVista guarantees the quality of its products as long as they are used according to their respective instructions and in accordance with national and international references. bioBoaVista does not take responsibility for the use of its products for purposes other than those described and approved by the company. All clinical diagnoses should be analyzed in conjunction with clinical evidence and not solely based on laboratory results.